




## Rapid & Reliable while maintaining the Versatility & Performance of a Traditional Wright-Giemsa Stain

QUICK III is a rapid hematology staining system utilizing three prepared solutions. The essential eosin and thiazine staining solutions are manufactured using certified dyes and quality control tested on human blood smears. Additionally, QUICK III has been formulated to provide superior results when used as an andrology stain. Users will be pleased with the differentiation of white blood cells and immature sperm with excellent resolution achieved in the staining of the head, neck and tailpiece.

- Quick 15 second slide dip procedure vs other products that may take up to 2 min.
- Reproducible and reliable results
- Blood smear results comparable to Wright-Giemsa protocols
- Provides superior results when used as an andrology stain
- Meets the specific needs of laboratories by offering 4oz, 16oz, and 1 gallon sizes



Staining Results	ASTRAL DIAGNOSTICS Quick III	Competitor's- 3 step Differential Stain
Red blood cells	Good staining	Good staining
Neutrophil	Excellent nuclear and cytoplasmic definition	Weak nuclear and cytoplasmic staining, lack of granular detail
Eosinophil	Excellent staining intensity; great nuclear detail	Weak staining, poor nuclear detail
Monocyte/Lymphocyte	Great staining detail	Good staining detail
Basophil	Great intensity staining	Inconsistent staining
Platelet	Good staining	Good staining
Immature sperm	Good staining	Good staining



Cat. #	Description	Size
5303	Quick III™ Stat-Pak	3 x 4oz
5316	Quick III Set	3 x 16oz
5317-16	Quick III Fixative	16oz
5317-G	Quick III Fixative	Gallon
5318-16	Quick III Solution I	16oz
5318-G	Quick III Solution I	Gallon
5319-16	Quick III Solution II	16oz
5319-G	Quick III Solution II	Gallon

Order today at [AstralDiagnostics.com](http://AstralDiagnostics.com)